

Clean Ears



Clean Ears

NAVEH PHARMA

1 Ver.3 04/20201

15mL Medical Device

For Ear Wax Removal Wash & Dissolve Jet Spray

- For ear wax removal
- Double action Wash & Dissolve system
- Prevents ear plugging
- Safe & effective
- User friendly
- Clinically tested
- No preservatives

Warning:

- CleanEars should not be used in cases of perforated eardrum and grommet
- Do not spray into the eye
- Keep out of reach of children



FOR EXTERNAL USE ONLY



Clean Ears

NAVEH PHARMA

15mL

For Ear Wax Removal Wash & Dissolve Jet Spray

Ingredients:

Mineral Oil, Phytosqualan, Spearmint Oil.



Directions:

Spray 1-2 puffs in each ear 3 times a day until the wax dissolves.

Storage Conditions:

Store in a cool place, do not expose to temperature over 50°C.



Naveh-Pharma LTD.

19, Yad-Harutzim St., Netanya, 4250519 Israel www.navehpharma.com



OBELIS S.A.

53 Boulevard Général Wahis 1030 Brussels, Belgium

www.navehpharma.com

LOT



7 290010 207091