



# Nasal Moisturizing Gel for Dry Nose

## An Advanced Moisturizing Gel for Dry Nasal Passages

**SalinaGel®** is an isotonic water-based saline nasal moisturizing gel that relieves dry and irritated nasal passages caused by aging, low humidity as well as side effects from medications for allergies, colds and acne.

Under certain conditions the ability of the nasal epithelial cells to secrete mucus, preserve moisture and filter inhaled air, may be impaired. Conditions such as low humidity, aging and side effects of medications to treat allergies, common colds and acne result in dryness in the nose and accumulation of crust causing discomfort, pain or coughing.

### **SalinaGel® Action**

**SalinaGel®** preserves nasal moisture and prevents dry and irritated nasal passage-ways caused by the daily exposure to changing air conditions, allergies and colds.

**SalinaGel®** achieves the following triple action:

- Moisturizes nasal airways and softens crust.
- Relieves dryness and irritation.
- Restores the filter and humidifying functions of the mucosa cells.

**SalinaGel®** ensures the normal function and wellness of the nasal cavities.

**SalinaGel®** is composed of sea-salt in isotonic concentration (similar to the salt concentration of the healthy tissue) with advanced moisturizing agents that preserve the healthy function of mucosa and restores the system's ability to filter inhaled air through the nose.

### **SalinaGel® Advantages**

- Consists of an advanced moisturising complex of Hyaluronic acid.
- Isotonic balance avoids discomfort and burning sensation.

### **Directions**

Gently blow your nose. Apply a small amount of gel directly with the nozzle or finger into the nostril and rub gently. Apply as needed.



**NAVEH**  
PHARMA

[www.navehpharma.com](http://www.navehpharma.com)